

Atomizer Nebulizer - Short User Guide

This PDF is a practical web version prepared for the landing page. For medical use, exact procedures and safety instructions, always follow the official manufacturer manual and the recommendation of a healthcare professional.

Safety note

- Use the device according to the official instructions and advice from a doctor or pharmacist.
- Do not add solutions or medication without advice from a qualified professional.
- Place the device on a stable, flat and dry surface.
- After use, clean and fully dry the accessories that come into contact with the solution.

Basic steps

- 1 Wash your hands and prepare a clean working surface.
- 2 Connect the air tube to the compressor and the nebulizer cup.
- 3 Add the recommended quantity of solution as advised by a professional.
- 4 Choose the adult mask, child mask or mouthpiece.
- 5 Turn on the device and use it according to the official instructions.
- 6 When finished, switch the device off and clean the accessories.

Technical data

Type	Piston compressor
Maximum pressure	185 kPa
Operating pressure	45 / 60 kPa
Air flow	8 / 11 l/min
Noise level	56 dBA
Nebulizer cup capacity	12 ml
Weight	1.4 kg
Warranty	3 years

Note: This document is a web guide. For exact medical procedures, the official manufacturer documentation applies.